


 Load MRHead  62-Probes  Start Navigation  Fast E-Field

[Phase 1] Surface Registration

- [Agent] Target loaded. Now, use the Stylus (or Mouse Click in simulation) to capture the **Nasion (Nose)**.
- [Captured] Nasion at -167.4, -109.1, -24.7
- [Agent] Good. Now capture the **Left Outer Canthus**.
- [Captured] Left Eye at -167.4, -109.1, -24.7
- [Agent] Good. Now capture the **Right Outer Canthus**.
- [Captured] Right Eye at -167.4, -109.1, -24.7
- [Result] Surface Registration Complete (Error: 62.36mm)



> Prompt: "Move coil to FC4 on the skin_surface"

> Expected Action: Agent automatically recognizes FC4 node, auto-corrects, and executes SlicerMCP pipeline.